# REDMOND

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**EN:** For your own safety read these instructions carefully before using the appliance.

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Complete set.....	4
Product specification.....	4
Model arrangement.....	6
I. BEFORE OPERATING.....	8
Safety measures.....	8
Prior to the first use.....	9
II. OPERATING OF A SCALES.....	9
Simple weighing mode.....	9
Weighing mode with constitution parameters indication.....	11
III. CLEANING AND MAINTENANCE.....	13
IV. BEFORE CONSULTING THE SERVICE CENTER.....	13
V. GUARANTEE.....	14

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**COMPLETE SET**

- |                      |       |
|----------------------|-------|
| 1. Electronic scales | 1 pc. |
| 2. Instruction       | 1 pc. |
| 3. Service booklet   | 1 pc. |

**TECHNICAL DATA**

Model	RS-713
Appliance type	electronic scales with constitution parameters
Power element	CR2032
Measurement range	1 - 150 kg
Measurement interval	50 g
Number of highly sensitive sensors	4
Electronic sensor control	+
Display illumination	+
Overload indicator	+
Unit of measure selection	+
Memory cells quantity	8
Automatic turn-on function	+
Automatic turn-off function	+

**THE MANUFACTURER HAS THE RIGHT OF MAKING MODIFICATIONS TO DESIGN, DELIVERED SET AND TECHNICAL DATA OF THE DEVICE IN THE COURSE OF IMPROVEMENT OF THE PRODUCTS WITHOUT ADDITIONAL NOTIFYING OF SUCH MODIFICATIONS.**

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
1. Scales platform
2. Number of user memory cell Indication
3. Age of the user input indication
4. Analytical mode input indication
5. Weight indication
6. Indication of a unit of measure of growth at input of parameters of the user
7. Indication of a unit of measure of weight
8. FAT mode indication – maintenance of a fat tissue
9. Indication of parameters of a constitution
10. TBW mode indication – Moisture indication
11. MUS mode indication – muscular tissue indication
12. Indication of a condition of a constitution (Thin, normal, plump, very plump)
13. Sensor button ▼ (DOWN) – for parameters selection
14. Sensor button SET – for entered parameters set
15. Sensor button ▲ (UP) – for parameters selection
16. Индикация режима BONE – bones tissue indication

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*Dear Customer!*

*Thank you for preferring the domestic electric appliance manufactured by **REDMOND**.*

***REDMOND** means quality, reliability and thoughtful attitude towards our clients.*

*We hope that you will select products of our firm in future, and the brand **REDMOND** will become predominant among your domestic appliances.*

***REDMOND** scales are high high-accuracy appliances for those leading healthy way of life and watching their weight carefully.*

## I. BEFORE OPERATING

### SAFETY MEASURES

**WHILE OPERATING THE APPLIANCE, FOLLOW THE INSTRUCTIONS GIVEN BELOW. IF USED PROPERLY, IT WILL SERVE YOU FOR MANY YEARS.**

- The appliance is intended for domestic use only.  
**The appliance is not intended for industrial use.**
- Never let children to use the appliance until they have been briefed and got the knowledge about the instructions of safe operation.  
**The appliance itself and the packing materials are not allowed to be used as children's playthings.**
- It is prohibited to perform independent repairs or make modifications to the design of the appliance.  
**All the service and repair works should be performed by authorised service-centre. Amateurish works performed can result in breaking, injuries and damage to property.**
- Use batteries of appropriate type only.  
**Battery type can be found in technical specifications or in the nameplate available on the appliance.**
- It is prohibited to immerse the appliance in water or any liquids.  
**This may result in damages.**

**OPERATION OF THE SCALES SHALL BE PROHIBITED IN THE FOLLOWING CASES**

- In premises with increased humidity, for example in a bath-room.  
**There is a risk of failure of the appliance.**
- Outdoors.  
**A risk exists, of malfunction due to rain and other weather conditions.**
- A mobile phone may create interference, which may lead to inaccurate results of weighing.
- Keep turned-on mobile phone at a distance of at least 2 metres from the scales or turn the mobile phone off.

**IT IS PROHIBITED TO USE THE APPLIANCE IN CASE OF ANY MALFUNCTIONS. THE MANUFACTURER SHALL NOT BE RESPONSIBLE FOR DAMAGE RESULTING FROM NON-OBSERVANCE OF THE REQUIREMENTS OF ACCIDENT PREVENTION AND THE RULES OF OPERATION OF THE ARTICLE.**

**KEEP THIS INSTRUCTION FOR USE IN THE FUTURE**

## **PRIOR TO THE FIRST USE**

1. Take out the appliance from the box carefully; remove all packing materials including advertising stickers.
2. Wipe the case of the appliance with moist cloth and then dry it up.
3. Remove the lid of the battery compartment and make sure the batteries are inserted correctly, with observance of polarity.
4. If the batteries have not been inserted, insert them in place.

**WHEN INSERTING THE BATTERIES, THE HANDS AND THE APPLIANCE SHOULD BE DRY.**

Connect the batteries to the connectors, observing the polarity.

5. Having connected the batteries, close the compartment lid.


## **II. OPERATING OF A SCALES**

### **SIMPLE WEIGHING MODE**

1. Select the desirable unit of measurement, by depressing the button kg/st/lb:  
**st** – stones  
**lb** – pounds  
**kg** – kilograms
2. Place the scales onto a level, solid horizontal surface.  
This condition is obligatory – it allows to receive exact results.

3. Scales is supplied with automatically switched on system, to switch on the scales step on scales platform.

**NOTE!**

*To get an exact value of weight it is necessary to observe a number of conditions:*

1. *Place the scales on flat hard horizontal surface.*
2. *Get on platform of scales smoothly.*
3. *Do not transfer the centre of gravity from one foot on another.*
4. *Avoid sharp movements and rockings during weighing.*

4. On LED display value of your weight will be displayed.

5. On the weighing termination – descend from a platform of scales. Scales will be disconnected automatically.
6. On LED display there will be a message  at in overload (over 150 kg)

### **WEIGHING MODE WITH CONSTITUTION PARAMETERS INDICATION**

7. Select the desirable unit of measurement, by depressing the button kg/st/lb:
  - st** – stones
  - lb** – pounds
  - kg** – kilograms
8. Place the scales onto a level, solid horizontal surface.  
This condition is obligatory – it allows to receive exact results.

**NOTE!**

*To get an exact value of weight it is necessary to observe a number of conditions:*

1. *Place the scales on flat hard horizontal surface.*
2. *Get on platform of scales smoothly.*
3. *Do not transfer the centre of gravity from one foot on another.*
4. *Avoid sharp movements and rockings during weighing.*

9. Press **SET** button, after that use SET and ▼ buttons, set number of memory cell (from 1 to 8), in which want, that your personal data remained.

**P1** – cell N1

**P2** – cell N2

.....

**P8** – cell N8

Press **SET** to confirm input and convert to next parameter.

4. Select analytic mode using ▲ and ▼ buttons:

- men mode


- women mode

- men mode who are regularly going in for sports


- women mode who are regularly going in for sports


Press **SET** to confirm input and convert to next parameter.

5. Enter Your height in cm using ▲ and ▼ buttons.


Press **SET** to confirm input and convert to next parameter.

6. Enter Your age using ▲ and ▼ buttons:


Press **SET** to confirm input and convert to next parameter.

7. After Your all personal data input, wait, while there will be zero on display.

8. Step on the scales – on the display your weight will be displayed.


9. After on the display value of your weight will be fixed, analytical process by definition of parameters of your constitution will begin.

10. Upon termination of analytical process, on the display parameters of your constitution will be displayed:

11. Upon termination of display of parameters of a constitution scales will be disconnected automatically.

### III. CLEANING AND MAINTENANCE

1. Never use abrasive cleaning agents, metallic brushes, spirit-based cleaning agents, etc. for clearing of the appliance.
2. It is prohibited to immerse the appliance into water, or wash under a flush of water.
3. Clear the platform of the scales with a damp and soft fabric, then wipe dry.

### IV. BEFORE CONSULTING THE SERVICE CENTER

Trouble	Probable cause	Remedy
The appliance would not operate.	The appliance has not been switched ON.	Switch ON the appliance, having stepped onto the platform of the scales. Replace power elements according to chapter «PRIOR TO THE FIRST USE».

Display readings do not correspond to the true value of weight.	Incorrect unit has been selected for measurement of weight.	Choose the necessary unit of weight measure using switcher. To get an exact value of weight it is necessary to observe a number of conditions: 1. Place the scales on flat hard horizontal surface. 2. Get on platform of scales smoothly. 3. Do not transfer the centre of gravity from one foot on another. 4. Avoid sharp movements and rockings during weighing.
The appliance turns off automatically.	Automatic turning-off system is activated.	This is normal. The scales are turned off automatically upon completion of operation of the same.

## V. GUARANTEE

The period of guarantee for the product is 1 year from the selling date.

The guarantee covers defects of material and workmanship.

Under the present guarantee, the manufacturer undertakes to repair or replace any part recognized as faulty on condition that the product is sent to the seller.

The present guarantee shall be valid only in case the product was used in compliance with this instruction, no modifications were made to it, it was not repaired or disassembled by unauthorized persons, or was not damaged due to improper handling. The complete set of the product must be retained also.

This guarantee does not cover natural deterioration of the product as well as brittle objects such as ceramic components, lamps, etc.

The guarantee shall not become effective unless the date of purchase is confirmed by the seller's signature and seal on the guarantee coupon.

Any repair under this guarantee shall not extend the period of guarantee and shall not initiate the beginning of a new guarantee.

**In compliance with paragraph 2, article 5 of the Federal Law of the RF «On protection of consumers' rights», the manufacturer has established a minimum service life for this product, which is 2 years from the date of manufacture on condition that the product is used in strict compliance with this manual and technical requirements.**

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